



# Avert<sup>®</sup> Sporocidal Disinfectant Cleaner

A ready-to-use disinfectant cleaner that also kills *Clostridium difficile* (*C. diff*) spores on hard, non-porous surfaces such as those found in hospitals and other healthcare facilities.

## Features & Benefits

- Effective in 4 minutes against *Clostridium difficile* (*C. diff*) spores. Cleans, disinfects, and deodorizes in one labor-saving step
- Bleach-based formula means no quat-binding concerns. 1 minute contact time effective against a broad spectrum of pathogens including bacteria, antibiotic resistant bacteria, viruses, and fungi
- Meets minimum active ingredient requirements for CDC, SHEA, APIC and AORN

## Applications

- Hard surface cleaning and disinfection
  - Designed for disinfecting medical hard, nonporous inanimate surfaces including: glass, glazed porcelain and stainless steel. Not recommended for use on: aluminum, brass, chipped enamel, clear plastic, copper, finished floors, natural marble, painted surfaces, and silver.
- Note: This product is compatible with the listed materials. If product is intended to be used on other surfaces, it is recommended to test product by applying to a small, discrete area to determine compatibility before proceeding with its use





# Avert<sup>®</sup> Sporicidal Disinfectant Cleaner

## Use instructions

### Special Instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores:

#### Personal Protection:

- Wear appropriate barrier protection such as gloves, gowns, masks or eye covering.

#### Cleaning Procedure:

- Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with the sporicidal product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.

#### Infectious Materials Disposal:

- Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

#### Disinfection Procedure:

- Apply product by coarse trigger sprayer to hard, nonporous inanimate surface.
- Spray 6-8 inches from surface, making sure to wet surfaces thoroughly.
- Surfaces must remain wet for 4 minutes.
- Air Dry, wipe surfaces to dry and remove any residue, or rinse with potable water as necessary.

#### For Use as a One-Step Cleaner/Disinfectant, Deodorizer, and Fungicide:

- Pre-clean heavily soiled areas.
- Apply product by coarse trigger sprayer to hard, nonporous inanimate surface.
- Spray 6-8 inches from surface, making sure to wet surfaces thoroughly.
- Wipe or rub with cloth, disposable wiper, sponge, brush, or hand mop to spread product on surface.
- Surfaces must remain wet for 1 minutes.
- Air Dry, wipe surfaces to dry and remove any residue, or rinse with potable water as necessary.

Technical data	Avert <sup>®</sup> Sporicidal Disinfectant Cleaner
EPA Reg. No.	70627-72
Color/Form	Pale Yellow, Liquid
pH	12.5
Scent	Chlorine
Shelf Life	12 months

Product	Pack size	Dilution	Product code
Avert <sup>®</sup> Sporicidal Disinfectant Cleaner	12 x 32 oz. / 946 mL Spray Bottles	Ready-to-use	100842725 

## Safe handling

Please make sure your employees read and understand the product label and Safety Data Sheet before using this product. The label contains directions for use; and both the label and SDS contain hazard warnings, precautionary statements and first aid procedures. SDS are available online at [www.diversey.com](http://www.diversey.com) or by calling 888.352.2249. Improper use or dilution may result in damage to surfaces and may result in health and physical hazards that match those of the concentrate. Please refer to the Diversey HazMat Library, only available through Internet Explorer, <http://naextranet.diversey.com/dot/>, for up-to-date shipping information.